


LAMP & Colormetric based COVID-19 detection kit

Rapid and sensitive visual detection of severe acute respiratory syndrome coronavirus 2 (SARS-COV-2) by reverse-transcription loop-mediated isothermal amplification with calcein and hydroxynaphthol blue


• Loop-mediated isothermal amplification (LAMP)



Loop-mediated isothermal amplification (LAMP) uses 4-6 primers that recognize 6-8 regions of target DNA.

DNA polymerase, which displaces the strand, starts synthesis and two of the primers form a loop structure to facilitate subsequent rounds of amplification.

• Colormetric : Calcein + Hydroxynaphthol blue



-Calcein-

pH 12 ↑ : Formation of a yellow-green chelation with Mg^{2+}
 Combined with Mn^{2+} : Black color(Quenching)
 Combined with Mg^{2+} : Yellow-green color


-HNB-

pH 12~13 : Changes from purple to blue
 when Mg^{2+} is present
 Combined with Mg^{2+} : Purple color
 Non-combined with Mg^{2+} : Blue color

LAMP & Colormetric based COVID-19 detection kit

Rapid and sensitive visual detection of severe acute respiratory syndrome coronavirus 2 (SARS-COV-2) by reverse-transcription loop-mediated isothermal amplification with calcein and hydroxynaphthol blue

- PCR results




Copies/5ul	Results
	Ct Mean
1×10^7	23.06
1×10^6	24.24
1×10^5	25.46
1×10^4	27.78
1×10^3	30.92
1×10^2	31.95
1×10^1	34.11
1×10^0	45.37
PC	22.01
NC	N/A

As a result of the PCR sensitivity test, it can be detected even at the level of 1 copy

LAMP & Colormetric based COVID-19 detection kit

Rapid and sensitive visual detection of severe acute respiratory syndrome coronavirus 2 (SARS-COV-2) by reverse-transcription loop-mediated isothermal amplification with calcein and hydroxynaphthol blue

- Colormetric results



Distinction between positive and negative colors
&
clear color change even in low-concentration samples